## **INKUBATOR 1000**

**new generation**

**Betriebsanleitung  
Instruction Manual  
Mode d'Emploi  
Instrucciones de  
Empleo  
Istruzioni per l'uso**

D	DEUTSCH	Seite 3 – 14
E	ENGLISH	page 15 – 26
F	FRANCAISE	page 27 – 38
ES	ESPAÑOL	página 39 – 50
I	ITALIANO	Pagina 51 – 62

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#### Niederspannungs-Richtlinie (73/23/EWG):

EN 61010-1 + EN 61010-2-010

Thank you for purchasing a Heidolph Instruments product. This item has been designed, made and inspected in compliance with DIN EN ISO 61010 for long-term and flawless operation.

## SUMMARY

<b>SUMMARY</b> .....	15
<b>STANDARD HARDWARE &amp; OPTIONS</b> .....	16
<b>GENERAL INFORMATION</b> .....	16
<b>SAFETY INFORMATION</b> .....	17
<b>SET-UP</b> .....	17
<b>OPERATION AND CONTROLS</b> .....	21
<b>CLEANING &amp; SERVICING</b> .....	22
<b>DISASSEMBLY &amp; STORAGE</b> .....	23
<b>DISPOSAL</b> .....	23
<b>TROUBLESHOOTING</b> .....	23
<b>SPECIFICATIONS</b> .....	24
<b>WARRANTY, LIABILITY &amp; COPYRIGHT</b> .....	25
<b>QUESTIONS / REPAIR WORK</b> .....	25
<b>CE-DECLARATION OF CONFORMITY</b> .....	26

**Important information**

**Advice about power cord / mains supply**

**Caution: mandatory action**

**Caution: fire- and explosion hazard**

**Advice about maintenance / repair**

## STANDARD HARDWARE & OPTIONS

	product	quantity	P/N 230/240V 50/60Hz	P/N 115V 50/60Hz
--	---------	----------	-------------------------	---------------------

### INCUBATOR 1000

	heater module	1	549-90010-00	549-90020-00
with	incubator hood, flat, consisting of: upper half, lower half, deflector plate, flat	1	549-90040-00	549-90040-00
or	incubator hood, high, consisting of: upper half, lower half, deflector plate, high	optional	549-90030-00	549-90030-00
	holder	1	21-005-014-48	21-005-014-48
	screw (to secure holder)	2	02-01-04-02-24	02-01-04-02-24
	spring washer (to secure holder)	2	02-03-01-01-09	02-03-01-01-09
	knurled bolt (to secure lower half of hood)	3	11-008-003-99	11-008-003-99
	Deflector plate	1	23-08-02-04-04	23-08-02-04-04
or	Deflector plate	1	23-08-02-04-03	23-08-02-04-03
	Instruction Manual	1	01-005-004-23	01-005-004-23
	Power cord	1	14-007-003-81	14-007-003-89

## GENERAL INFORMATION

Unpack your item carefully.

Inspect for damage and report such damage or missing parts to your supplier right away.

Read your Instruction Manual carefully. Take time to save time when working with your product. Make sure that every user has read and understood the Instruction Manual.

Please store the Instruction Manual in a place easily accessible to every user.

**IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !**

A so-called EURO-plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V) is standard on all of the products.

For the Continental US they feature a US-standard plug (NEMA Pub.No.WDI.1961 ASA C 73.1 . 1961 page 8 15A 125V).

For using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one suiting your needs.

As shipped, the item features a protective ground wire. When replacing the original plug, make sure to reconnect this protective ground wire in the new plug !

## SAFETY INFORMATION

**Please comply with all safety and accident-prevention regulations as in force for laboratory work!**

**Use extra care when working with flammable substances; refer to safety data sheets.**

**Use extra care when working in the vicinity of flammable and explosive substances. Motors are of non-sparking type, the item itself however is not explosion-proof.**

**When connecting your item with your local power supply, please make sure your item is designed for your local voltage; go by the data plate on the item.**

**Please connect your item with a protective-ground outlet only.**

**Turn your power switch OFF whenever the item is not used, or before disconnecting the plug.**

**Repair work is limited to technicians approved by Heidolph Instruments.**

**Your item needs a solid stand.**

**Lab bench needs to be of rigid design, and features an anti-skid surface coat.**

**Before starting the item, make sure all samples are attached safely (must not move while shaking).**

## SET-UP

The INCUBATOR 1000, consists of heater module and hood; it may be combined with shakers and mixers as stipulated below. The incubator heats fluids in whatever sort of vessel during shaking. The hood is made of clear plastic material to allow watching of vessels during shaking.

**The INCUBATOR 1000 may be combined with shakers and mixers stipulated below:**

TITRAMAX 1000

UNIMAX 1010

PROMAX 1020


DUOMAX 1030

POLYMAX 1040

Preparing shakers and mixers for use with / without the incubator can be done in a few steps and will take only a few minutes.

# E


- ♦ Install holder (E) with two screws and spring washer from your hardware bag on bottom of heater module (B).

- ♦ Locate heater module on the right of shaker (A).


- ♦ Plug rubber feet of shaker in cutouts of holder.
- ♦ Install lower part of hood (C) in shaker and heater module.  
If necessary move tray so that cut outs of the lower part of hood (C) can be joined with the tray.


**E**

- ◆ The third attachment point is placed on the left hand side of the apparatus (N). Turn third knurled screw (G1) into thread (N) on left hand side and tighten. Lower part (C) is now secure.

# E

- ◆ Secure lower half (C) on heater module; use 2 ea. knurled bolts from your hardware bag.

- ◆ Place hood (D) on lower part. To do so loosen fastening screw (H) by turning counter-clockwise until gap is ca. 6mm. Align hinges to an upright position. While mounting hood (D) place threaded part of fastening screw (H) in slash of hinges (J).
- ◆ Tighten fastening screws (H).
- ◆ Open hood and plug deflector plate (F) in cutouts of lower half (C). This plate is important for proper distribution of hot air under the hood.
- ◆ heater module power supply is independent of shaker and mixer. Heating module comes with its own power cord. Use handle (L) to open/close hood (D).
- ◆ When closed, the integrated magnets (M) keep the hood (D) in position.

- The incubator has friction hinges. It allows the user to open it to any position. System related, a moderate force is required to open and close the hood.
- For cleaning, hood can be dismounted. To do so undo fastening screws (H) and lift off hood.

## OPERATION AND CONTROLS

- Before connecting power cord with mains outlet, make sure that:
  - your item is designed for your local voltage and frequency (data plate on item).
  - master switch is set to "0" and all controls are in "min." position (turn CCW completely); this way you avoid spilling fluids by too intense shaking action.
- Carefully close your vessels and select appropriate shaking intensity (if vessel remain open) to avoid splashes and spillage.
- We recommend to start with low shaking speed and gradually increase frequency to avoid accidental shaking at high speed.
- When handling hazardous fluids, make reference to applicable safety information.
- Do not close ventholes of heater module.  
Maintain unobstructed circulation of air under hood..
- Avoid entry of foreign matter or spilled fluids through ventholes in heater module.  
Negligence might damage fan (heater element).

### Incubator 1000, heating module

The front panel features the following controls (from right to left):

# E

- A 2-pole master switch, lighted green
- B temperature control
- C heater LED
- D LED-display: actual temperature
- E LED-display: nominal temperature

- ◆ close master switch to start the incubator. Fan runs continuously to circulate air under the hood.
- ◆ adjust nominal temperature with temperature control (rotating knob); this value is displayed continuously by the lefthand, 2-digit LED display (E).
- ◆ actual temperature is displayed continuously by the righthand, 2-digit LED display (D).
- ◆ yellow LED (C) is lighting when item is heating.
- ◆ the heater module features overheat protection. Sensor is located on the heater element. In case of malfunction, heater turns off automatically. The fan continues running, all displays remain on. After temperature settling, the heater is turned on again.

Transfer of heat from circulating air to beakers or flasks might require some time; we recommend to pre-heat cold flasks or beakers.

Working with aqueous substances at elevated temperatures will involve a certain level of evaporation from open beakers or flasks.

Accuracy of incubator's temperature control will be reached at nominal temperatures of >5K above ambient only.

## CLEANING & SERVICING

**Cleaning:** wipe housing clean with a damp cloth (add some sort of mild liquid soap).

### Note

To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based detergents, abrasive substances, ammonia, rags or cleaning agents containing metal particles.

The item is maintenance-free. Repair work is limited to technicians approved or appointed by Heidolph Instruments. Please call your local Heidolph Instruments Dealer or a Heidolph Instruments Field Representative (also refer to page 26)

## DISASSEMBLY & STORAGE

### **Disassembly**

1. Turn item OFF and disconnect mains plug.
2. Remove all of the hardware arranged around the incubator to ease disassembly.
3. Unload all flasks and beakers from incubator, uninstall optional equipment.

### **Forward & Store**

1. We recommend to store the item and its components in its original box, or a similar container that offers adequate protection against damage in transit. Tape the box securely.
2. Store the item in a dry place.

Caution

Do not jolt or shake the item during transport.

## DISPOSAL

For disposal, please comply with your local or national regulations.

Split by metal, plastic, etc.

Packing material to be treated as described above (material split).

## TROUBLESHOOTING

Work on electric and electronic components is limited to qualified personnel.

### **Heater module master switch dim**

1. check power cord
2. check circuit breakers on back side

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### **Heater module won't heat**

1. nominal temperature set to a value below actual ambient.
2. overheat protection disconnected heater.  
Remedy:  
allow heater to cool down, check vent holes for clogging.
3. electronic problem (annunciators dim).

## SPECIFICATIONS

### **INCUBATOR 1000**

power supply	230/240 V 50/60 Hz / 115 V 50/60 Hz (optional)
temperature range	5 K above ambient (65°C max.)
heating power	300 W
accuracy (temperature control)	± 2 K up to 50°C ± 4 K above 50°C (air temperature)
temperature distribution	± 2 K up to 50°C ± 3 K above 50°C
time to heat circulating air	abt. 12 minutes from ambient to 50°C abt. 25 minutes from ambient to 60°C (unobstructed circulation)
protection class	IP 30
duty cycle	100 %
temperature control	electronic
temperature sensor	Pt 100, installed in heater module
heater	tubular heater with overheat protection
circulation of air	cross current fan, recirculating principle
temperature display	nominal: 2-digit LED display actual: 2-digit LED display
size (w x d x h)	abt. 170 mm x 375 mm x 97 mm

### **Hood for Incubator 1000**

material	PETG, split-design
size (w x d x h)	abt. 570 mm x 390 mm x 270 mm
volume (air)	abt. 50 dm <sup>3</sup>

### **Hood (flat) for Incubator 1000**

material	PETG, split-design
size (w x d x h)	abt. 570 mm x 390 mm x 170 mm
volume (air)	abt. 30 dm <sup>3</sup>

## WARRANTY, LIABILITY & COPYRIGHT

### **Warranty**

Heidolph Instruments warrants that the present product shall be free from defects in material (except wear parts) and workmanship for 3 years from the date shipped off the manufacturer's warehouse.

Transit damage is excluded from this warranty.

To file for such warranty service, contact Heidolph Instruments (phone:++49-9122-9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty.

Alterations to the present warranty need Heidolph Instruments' consent in writing.

### **Exclusion Clause**

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

### **Copyright**

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

## QUESTIONS / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please contact the following address:

For repair services please call Heidolph Instruments (phone: +49 - 9122 - 9920-68) or your local, authorized Heidolph Instruments Dealer.

### **Note**

You will receive approval for sending your defective item to the following address:

**Heidolph Instruments GmbH & Co. KG**  
**Lab Equipment Sales**  
**Walpersdorfer Str. 12**  
**D-91126 Schwabach / Germany**  
**Tel.: +49 – 9122 - 9920-68**  
**Fax: +49 – 9122 - 9920-65**  
**E-Mail: Sales@Heidolph.de**

### **Note**

If you are based in the United States of America, please contact Heidolph US:

**Heidolph Instruments, LLC**  
**Lab Equipment Sales**  
**2615 River Rd.**  
**Cinnaminson, NJ 08077**  
**Phone: 856-829-6160**  
**Fax: 856-829-7639**  
**E-Mail: heidolph@snip.net**

### Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
- take proper protective measure to ensure the safety of our receiving and service personnel
- mark the pack IAW Hazardous Materials Act

### CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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